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
Reference

Figure A — Important Components of an HPLC System

Isolating HPLC Problems

In an HPLC system, problems can arise from many sources. First define the problem, then isolate the source.

Use Table 1 to determine which component(s) may be causing the trouble. Thereafter, a process of elimination will usually enable you to pinpoint the specific cause and correct the problem.

How to Prevent Mobile Phase Problems

Low sensitivity and rising baselines, noise, or spikes on the chromatogram can often be attributed to the mobile phase. Contaminants in the mobile phase are especially troublesome in gradient elution. The baseline may rise, and spurious peaks can appear as the level of the contaminated component increases.

Water is the most common source of contamination in reversed phase analyses. You should use only high purity distilled or deionized water when formulating mobile phases. However, several common deionizers introduce organic contaminants into the water. To remove these contaminants, pass the deionized water through activated charcoal or a C18 cartridge.

Use only HPLC grade solvents, salts, ion pair reagents, and basic modifiers. Cleaning lower quality solvents is time consuming, and trace levels of contaminants often remain. These trace con-

taminants can cause problems when you use a high sensitivity ultraviolet or fluorescence detector.

To help prevent bubbles in the system, degas the mobile phase prior to use. Spurge it with helium (3-5psi) during use. Filtering the mobile phase through a 0.2 or 0.45 micron filter eliminates gas. This will also remove particles that could produce noisy baselines or plug the column.

Avoid mixing dissimilar or immiscible solvents and remember to flush the system when going from one mobile phase to another. To ensure complete flushing, always use a solvent of intermediate polarity (Table 3, page 15).

Because many aqueous buffers promote the growth of algae or bacteria, you should discard cloudy buffers and make them up fresh. Prevent microorganism growth by adding about 100ppm of sodium azide to aqueous buffers. Or you can mix these buffers with an organic solvent such as acetonitrile.

Use ion pair reagents carefully. The optimum chain length and concentration of the reagent must be determined for each analysis. In general, increasing the concentration or chain length increases retention times. We recommend using concentrations of 0.2 to 10mM. High concentrations (>50%) of acetonitrile and some other organic solvents can precipitate ion pair reagents.

Volatile basic modifiers, such as triethylamine, are useful when you wish

to recover a compound for further analysis. These modifiers also let you avoid problems associated with ion pair reagents. They can be added to the buffer at concentrations of 0.1 to 1.0%. Increasing the concentration improves peak shape for basic compounds, but can alter retention times.

Recent studies suggest that recycling mobile phase used for isocratic separations prolongs column life and reduces solvent costs. However, you must adjust the recorder or integrator to compensate for the gradual baseline increase that usually follows. Remember . . . volatile components can evaporate, especially when the mobile phase is recycled or continuously degassed, and this will change mobile phase composition.

Isolating Pump Problems

The pump must deliver a constant flow of solvent to the column over a wide range of conditions. Modern HPLC pump designs include syringe, diaphragm, piston, or combinations.

Pumping system problems are usually easy to spot and correct. Some of the more common symptoms are erratic retention times and a noisy baseline, or spikes on the chromatogram. Leaks at pump fittings or seals will result in poor chromatography. A sure sign of a leak is a buildup of salts at a pump connection.

To isolate and repair specific problems related to your apparatus, use the troubleshooting and maintenance sec-

tions of the operation manual. Pump seals may require periodic replacement. You should perform regular maintenance rather than waiting to deal with a problem.

Injector and Injection Solvents

The injector rapidly introduces the sample into the system with minimal disruption of the solvent flow. HPLC systems currently use variable loop, fixed loop, and syringe-type injectors. These are activated manually, pneumatically, or electrically.

Mechanical problems involving the injector (e.g., leaks, plugged capillary tubing, worn seals) are easy to spot and correct. Other problems, such as irreproducible injections, are more difficult to solve. But these, too, can sometimes be traced to the injector.

Variable peak heights, split peaks, and broad peaks can be caused by incompletely filled sample loops, incompatibility of the injection solvent with the mobile phase, or poor sample solubility. Whenever possible, dissolve and inject samples in mobile phase. Otherwise, be sure the injection solvent is of lower eluting strength than the mobile phase (Table 3).

Column Protection

Although not an integral part of most equipment, mobile phase filters, in-line filters, saturator columns, and guard columns greatly reduce problems associated with complex separations.

Filters and guard columns prevent particles and strongly retained compounds from accumulating on the analytical column. Silica particles in a saturator column dissolve in basic mobile phases, protecting the silica based packing in the analytical column.

The useful life of these disposable products depends on many factors, including mobile phase composition, sample purity, pH, etc. As these devices become contaminated or plugged with particles, pressure increases and peaks broaden or split. For more about column protection, see the products pages of this guide and request Bulletin 781.

Getting the Most from Your Analytical Column

Columns are available in many different sizes and designs. A wide variety of packings are also available, but all have the same purpose — to perform the separation. The most common problem associated with analytical columns

is deterioration. This is true regardless of whether the column contains a bonded reversed or normal phase, ion exchange, size exclusion, or resin based packing.

Symptoms of deterioration are poor peak shape, split peaks, shoulders, loss of resolution, decreased retention times, and high back pressure. These symptoms indicate contaminants have accumulated on the frit or column inlet, or there are voids, channels, or a depression in the packing bed.

Deterioration is more evident in higher efficiency columns. For example, a 3 micron packing retained by 0.5 micron frits is more susceptible to plugging than a 5 or 10 micron packing retained by 2 micron or larger frits. Proper column protection and sample preparation are essential to getting the most from each column.

Overloading a column can cause poor peak shapes and other problems. Column capacity depends on many factors, but typical values are:

Analytical column (25cm x 4.6mm)	< 500µg
Semi-preparative column (25cm x 10mm)	< 100mg
Preparative column (25cm x 21.2mm)	< 500mg

Solving Detector Problems

More than 20 types of detectors are available for HPLC systems. The most common are fixed and variable wavelength ultraviolet photometers and refractive index detectors. Electrochemical and conductivity detectors are currently popular. Improvements in detector cell technology have made them more durable and easier to use.

Detector problems fall into two categories — electrical and mechanical/optical. For electrical problems, you should deal with the instrument manufacturer. Mechanical or optical problems can usually be traced to the flow cell. Detector-related problems include leaks, air bubbles, and cell contamination. These usually produce spikes or baseline noise on the chromatograms.

Some cell materials — especially those used in refractive index detectors — are sensitive to high pressure. Flow rates or back pressures that exceed the manufacturer's recommendation will break the cell window. Old or defective lamps and incorrectly set detector rise time, gain, or attenuation will reduce sensitivity and peak height. Faulty or reversed cable connections can also cause problems.

Column Heater, Recorder

These components seldom cause problems with the system. They will be discussed in the troubleshooting table (Table 1).

Keep Accurate Records


Most problems don't occur overnight, but develop gradually. Accurate record keeping, then, is vital to detecting and solving many problems.



Evaluate every column you receive, when you receive it and at regular intervals thereafter. By keeping a written history of column efficiency, mobile phases used, lamp current, pump performance, etc., you can monitor your system's performance.



Records also help prevent mistakes, such as introducing water into a silica column or precipitating buffer in the system. Many analysts modify their HPLC systems in some way. Reliable records are the best way to ensure that a modification does not introduce problems. For problems relating to pumps, detectors, automatic samplers, and data systems, consult your instrument manual's troubleshooting guide.




PROBLEM INDEX					
Problem	Problem No.	Problem	Problem No.	Problem	Problem
Baseline —		Ghost peaks	19	Peaks —	
drift	12	Peak shapes, incorrect —		height change	16
noise, irregular	14	broad	15	missing	2
noise, regular	13	fronting	10	negative	18
Column back pressure —		rounded	11	no peaks	1
higher than usual	4	split	7	unresolved	6
lower than usual	3	tailing	8, 9	Retention times, variable	5
				Selectivity change	17




Table 1 — HPLC Problems, Probable Causes, and Remedies

Problem	Probable Cause	Remedy/Comments
Problem No. 1		
<p>No Peaks/Very-Small Peaks</p>  <p>Normal</p> <p>Problem</p> <p>Problem</p>	<ol style="list-style-type: none"> 1. Detector lamp off. 2. Loose/broken wire between detector and integrator or recorder. 3. No mobile phase flow. 4. No sample/deteriorated sample. 5. Settings too high on detector or recorder. 	<ol style="list-style-type: none"> 1. Turn lamp on. 2. Check electrical connections and cables. 3. See "No Flow" (Problem No. 2). 4. Be sure automatic sampler vials have sufficient liquid. Evaluate system performance with fresh standard to confirm sample as source of problem. 5. Check attenuation or gain settings.



Problem	Probable Cause	Remedy/Comments
<p>Problem No. 2</p> <p>No Flow</p>  <p>Normal</p> <p>Problem</p>	<ol style="list-style-type: none"> 1. Pump off. 2. Flow interrupted/obstructed/ 3. Leak. 4. Air trapped in pump head. (Revealed by pressure fluctuations.) 	<ol style="list-style-type: none"> 1. Start pump. 2. Check mobile phase level in reservoir(s). Check flow through-out system. Examine sample loop for obstruction or air lock. Make sure mobile phase inlet filter is clean. Make sure mobile phase components are miscible and mobile phase is properly degassed. 3. Check system for loose fittings. Check pump for leaks, salt build-up, unusual noises. Change pump seals if necessary. 4. Disconnect tubing at guard column (if present) or analytical column inlet. Check for flow. Purge pump at high flow rate (10ml/min.), prime system if necessary. (Prime each pump head separately.) If system has check valve, loosen valve to allow air to escape. If problem persists, flush system with 100% methanol. If problem still persists, contact system manufacturer.
<p>Problem No. 3</p> <p>No Pressure/Pressure Lower than Usual</p> 	<ol style="list-style-type: none"> 1. Leak. 2. Mobile phase flow interrupted/obstructed. 3. Air trapped in pump head. (Revealed by pressure fluctuations.) 4. Leak at column inlet end fitting. 5. Air trapped elsewhere in system. 6. Worn pump seal causing leaks around pump head. 	<ol style="list-style-type: none"> 1. Check system for loose fittings. Check pump for leaks, salt build-up, unusual noises. Change pump seals if necessary. 2. Check mobile phase level in reservoir(s). Check flow through-out system. Examine sample loop for obstruction or air lock. Make sure mobile phase inlet filter is clean. Make sure mobile phase components are miscible and mobile phase is properly degassed. 3. Disconnect tubing at guard column (if present) or analytical column inlet. Check for flow. Purge pump at high flow rate (10ml/min.), prime system if necessary. (Prime each pump head separately.) If system has check valve, loosen valve to allow air to escape. 4. Reconnect column and pump solvent at 1-2ml/min. If pressure is still low, check for leaks at inlet fitting or column end fitting. 5. Disconnect guard and analytical column and purge system. Reconnect column(s). If problem persists, flush system with 100% methanol. 6. Replace seal.

Problem	Probable Cause	Remedy/Comments
<p><i>Problem No. 4</i></p> <p>Pressure Higher than Usual</p> 	<ol style="list-style-type: none"> 1. Problem in pump, injector, or tubing. 2. Obstructed guard or analytical column. 	<ol style="list-style-type: none"> 1. Disconnect guard and analytical column, run pump at 2-5ml/min. If pressure is minimal, see Cause 2. If not, isolate cause by systematically eliminating system components, starting with detector and working back to pump. 2. Remove guard column (if present) and check pressure. Replace guard column if necessary. If analytical column is obstructed, reverse and flush (page 13). If problem persists, column may be clogged with strongly retained contaminants. Use appropriate restoration procedure (Table 2, page 14). If problem still persists, change inlet frit (page 13) or replace column.
<p><i>Problem No. 5</i></p> <p>Variable Retention Times</p> 	<ol style="list-style-type: none"> 1. Leak. 2. Change in mobile phase composition. (Small changes can lead to large changes in retention times.) 3. Air trapped in pump. (Retention times increase and decrease at random times.) 4. Column temperature fluctuations (especially evident in ion exchange systems). 5. Column overloading. (Retention times usually decrease as mass of solute injected on column exceeds column capacity.) 6. Sample solvent incompatible with mobile phase. 7. Column problem. (Not a common cause of erratic retention. As a column ages, retention times gradually decrease.) 	<ol style="list-style-type: none"> 1. Check system for loose fittings. Check pump for leaks, salt build-up, unusual noises. Change pump seals if necessary. 2. Check make-up of mobile phase. If mobile phase is machine mixed, hand mix and supply from one reservoir. 3. Purge air from pump head or check valves. Change pump seals if necessary. Be sure mobile phase is degassed. 4. Use reliable column oven or insulate column. (Note: higher column temperatures increase column efficiency. For optimum results, heat eluant before introducing it onto column.) 5. Inject smaller volume (e.g., 10μl vs. 100μl) or 1:10 and 1:100 dilutions of sample. 6. Adjust solvent. Whenever possible, inject samples in mobile phase. 7. Substitute new column of same type to confirm column as cause. Discard old column if restoration procedures fail (see page 14).

Problem	Probable Cause	Remedy/Comments
Problem No. 6		
<p>Loss of Resolution</p>  <p>Normal</p> <p>Problem</p>	<ol style="list-style-type: none"> 1. Mobile phase contaminated/deteriorated (causing retention times to change). 2. Obstructed guard or analytical column. 	<ol style="list-style-type: none"> 1. Check make-up of mobile phase. (page 2). 2. Remove guard column (if present) and attempt analysis. Replace guard column if necessary. If analytical column is obstructed, reverse and flush (page 13). If problem persists, column may be clogged with strongly retained contaminants. Use appropriate restoration procedure (Table 2, page 14). If problem still persists, change inlet frit (page 13) or replace column.
Problem No. 7		
<p>Split Peaks</p>  <p>Normal</p> <p>Problem</p>	<ol style="list-style-type: none"> 1. Contamination on guard or analytical column inlet. 2. Sample solvent incompatible with mobile phase. 	<ol style="list-style-type: none"> 1. Remove guard column (if present) and attempt analysis. Replace guard column if necessary. If analytical column is obstructed, reverse and flush (page 13). If problem persists, column may be clogged with strongly retained contaminants. Use appropriate restoration procedure (Table 2, page 14). If problem still persists, inlet frit is probably plugged. Change frit (page 13) or replace column. 2. Adjust solvent. Whenever possible, inject samples in mobile phase.
Problem No. 8		
<p>Peaks Tail on Initial and Later Injections</p>  <p>Normal</p> <p>Problem</p>	<ol style="list-style-type: none"> 1. Interference in sample. 2. Wrong column type. 3. Wrong mobile phase pH. 4. Sample reacting with active sites. 	<ol style="list-style-type: none"> 1. Check column performance with standards. 2. Try another column type (e.g., deactivated column for basic compounds). 3. Adjust pH. For basic compounds, lower pH usually provides more symmetric peaks. 4. Add ion pair reagent or volatile basic modifier (page 3).

Problem	Probable Cause	Remedy/Comments
<p><i>Problem No. 9</i></p> <p>Peaks Formerly Symmetric, Now Tailing</p> 	<ol style="list-style-type: none"> 1. Interference in sample. 2. Mobile phase contaminated/deteriorated. 3. Guard or analytical column contaminated/worn out. 	<ol style="list-style-type: none"> 1. Check column performance with standards. 2. Check make-up of mobile phase (page 2). 3. Remove guard column (if present) and attempt analysis. Replace guard column if necessary. If analytical column is source of problem, use appropriate restoration procedure (Table 2, page 14). If problem persists, replace column.
<p><i>Problem No. 10</i></p> <p>Fronting Peaks</p> 	<ol style="list-style-type: none"> 1. Interference in sample. 2. Shoulder or gradual baseline rise before a main peak may be another sample component. 3. Column overloaded. 4. Sample solvent incompatible with mobile phase. 	<ol style="list-style-type: none"> 1. Check column performance with standards. 2. Increase efficiency or change selectivity of system to improve resolution. Try another column type if necessary (e.g., switch from nonpolar C18 to polar cyano phase). 3. Inject smaller volume (e.g., 10 μl vs. 100 μl) or 1:10 and 1:100 dilutions of sample. 4. Adjust solvent. Whenever possible, inject samples in mobile phase. Flush polar bonded phase column with 200ml HPLC grade ethyl acetate (2.5ml/min.), then with intermediate polarity solvent prior to analysis.
<p><i>Problem No. 11</i></p> <p>Rounded Peaks</p> 	<ol style="list-style-type: none"> 1. Detector operating outside linear dynamic range. 2. Recorder gain set too low. 3. Column overloaded. 4. Sample-column interaction. 5. Detector and/or recorder time constants are set too high. 	<ol style="list-style-type: none"> 1. Reduce sample volume and/or concentration. 2. Adjust gain. 3. Inject smaller volume (e.g., 10 μl vs. 100 μl) or 1:10 and 1:100 dilutions of sample. 4. Change buffer strength, pH, or mobile phase composition. If necessary, raise column temperature or change column type. (Analysis of solute structure may help predict interaction.) 5. Reduce settings to lowest values or values at which no further improvements are seen.

Problem	Probable Cause	Remedy/Comments
<p>Problem No. 12</p> <p>Baseline Drift</p> <p>Normal</p> <p>Problem</p>	<ol style="list-style-type: none"> 1. Column temperature fluctuation. (Even small changes cause cyclic baseline rise and fall. Most often affects refractive index and conductivity detectors, UV detectors at high sensitivity or in indirect photometric mode.) 2. Nonhomogeneous mobile phase. (Drift usually to higher absorbance, rather than cyclic pattern from temperature fluctuation.) 3. Contaminant or air buildup in detector cell. 4. Plugged outlet line after detector. (High pressure cracks cell window, producing noisy baseline.) 5. Mobile phase mixing problem or change in flow rate. 6. Slow column equilibration, especially when changing mobile phase. 7. Mobile phase contaminated, deteriorated, or prepared from low quality materials. 8. Strongly retained materials in sample (high k') can elute as very broad peaks and appear to be a rising baseline. (Gradient analyses can aggravate problem.) 9. Mobile phase recycled but detector not adjusted. 10. Detector (UV) not set at absorbance maximum but at slope of curve. 	<ol style="list-style-type: none"> 1. Control column and mobile phase temperature, use heat exchanger before detector. 2. Use HPLC grade solvents, high purity salts, and additives. Degas mobile phase before use, sparge with helium during use. 3. Flush cell with methanol or other strong solvent. If necessary, clean cell with 1N HNO_3 (never with HCl). 4. Unplug or replace line. Refer to detector manual to replace window. 5. Correct composition/flow rate. To avoid problem, routinely monitor composition and flow rate. 6. Flush column with intermediate strength solvent, run 10-20 column volumes of new mobile phase through column before analysis. 7. Check make-up of mobile phase (page 2). 8. Use guard column. If necessary, flush column with strong solvent between injections or periodically during analysis. 9. Reset baseline. Use new materials when dynamic range of detector is exceeded (page 3). 10. Change wavelength to UV absorbance maximum.
<p>Problem No. 13</p> <p>Baseline Noise (regular)</p> <p>Normal</p> <p>Problem</p>	<ol style="list-style-type: none"> 1. Air in mobile phase, detector cell, or pump. 2. Leak. 3. Incomplete mobile phase mixing. 4. Temperature effect (column at high temperature, detector unheated). 5. Other electronic equipment on same line. 6. Pump pulsations. 	<ol style="list-style-type: none"> 1. Degas mobile phase. Flush system to remove air from detector cell or pump. 2. Check system for loose fittings. Check pump for leaks, salt buildup, unusual noises. Change pump seals if necessary. 3. Mix mobile phase by hand or use less viscous solvent. 4. Reduce differential or add heat exchanger. 5. Isolate LC, detector, recorder to determine if source of problem is external. Correct as necessary. 6. Incorporate pulse dampener into system.

Problem	Probable Cause	Remedy/Comments
<i>Problem No. 14</i>		
<p>Baseline Noise (irregular)</p> 	<ol style="list-style-type: none"> 1. Leak. 2. Mobile phase contaminated, deteriorated, or prepared from low quality materials. 3. Detector/recorder electronics. 4. Air trapped in system. 5. Air bubbles in detector. 6. Detector cell contaminated. (Even small amounts of contaminants can cause noise.) 7. Weak detector lamp. 8. Column leaking silica or packing material. 	<ol style="list-style-type: none"> 1. Check system for loose fittings. Check pump for leaks, salt build-up, unusual noises. Change pump seals if necessary. 2. Check make-up of mobile phase. (page 2). 3. Isolate detector and recorder electronically. Refer to instruction manual to correct problem. 4. Flush system with strong solvent. 5. Purge detector. Install back pressure device after detector. 6. Clean cell. 7. Replace lamp. 8. Replace column.
<i>Problem No. 15</i>		
<p>Broad Peaks</p> 	<ol style="list-style-type: none"> 1. Mobile phase composition change. 2. Mobile phase flow rate too low. 3. Leak (especially between column and detector). 4. Detector settings incorrect. 5. Extra-column effects: <ol style="list-style-type: none"> a. Column overloaded. b. Detector response time or cell volume too large. c. Tubing between column and detector too long or ID too large. d. Recorder response time too high. 6. Buffer concentration too low. 7. Guard column contaminated/worn out. 8. Column contaminated/worn out. 9. Void at column inlet. 10. Peak represents two or more poorly resolved compounds. 11. Column temperature too low. 	<ol style="list-style-type: none"> 1. Prepare new mobile phase. 2. Adjust flow rate. 3. Check system for loose fittings. Check pump for leaks, salt build-up, and unusual noises. Change pump seals if necessary. 4. Adjust settings. 5. <ol style="list-style-type: none"> a. Inject smaller volume (e.g., 10μl vs. 100μl) or 1:10 and 1:100 dilutions of sample. b. Reduce response time or use smaller cell. c. Use as short a piece of 0.007-0.010" ID tubing as practical. d. Reduce response time. 6. Increase concentration. 7. Replace guard column. 8. Replace column with new one of same type. If new column provides symmetrical peaks, flush old column (Table 2, page 14), then retest. 9. Open inlet end and fill void (page 13) or replace column. 10. Change column type to improve separation. 11. Increase temperature. Do not exceed 75°C unless higher temperatures are acceptable to column manufacturer.


Figure 1(a-c): Effect of cartridge length in isocratic elution under identical conditions.

Figure 2(a-c): Effect of cartridge length in gradient elution under identical conditions.

To Minimize Column Cost

The selection of the shortest cartridge that gives satisfactory resolution offers multiple advantages.

As shown in Figure 1, shorter cartridges reduce analysis time and increase sample detectability. Shorter cartridges are low in cost and also reduce solvent consumption. If you can resolve your sample on a short cartridge, you can save both time and money by replacing a long

cartridge that has more resolution than is required.

The effect of cartridge length on resolution is less in gradient elution than in isocratic separations, because solvent composition (not cartridge length) is the dominant force in gradient elution. Figure 2 shows that short cartridge columns produce resolutions nearly as good as long columns.

Low cost, easily interchanged, short Pierce cartridge columns (especially the 3 cm length) are ideal for methods development. You can afford to have a wide variety of column packing materials on hand for convenient, time-saving evaluations. And when the sorbent has been selected and higher resolution is required, simply choose a longer cartridge with the same packing.

Restoring Your Column's Performance

These procedures should rejuvenate a column whose performance has deteriorated.

Disconnect and reverse the column. Connect it to the pump, but not the detector. Follow the appropriate flushing procedure in this table, using a flow rate that results in column back pressure of 1500-4000psi. If you have a SUPELCOSIL® column, analyze the test mix provided using the conditions on the quality control test sheet. Efficiency, symmetry, and capacity should be within 10-15% of the test sheet values. If not, repack the column inlet (page 13) or replace the column.

Table 2 — Column Restoration Procedures

<p>Nonbonded Silica Column</p> <ul style="list-style-type: none"> Flush with 50ml hexane Flush with 50ml methylene chloride Flush with 50ml isopropanol Flush with 50ml methanol Flush with 25ml methylene chloride Flush with 25ml mobile phase Evaluate column performance 	<p>Polar Column (amino, cyano, or diol column, SUPELCOSIL® LC-PCN, or LC-(R)-Urea column)</p> <p>For a column used in the reversed phase mode (e.g., drug analyses using an organic solvent/aqueous buffer mobile phase), follow the same cleanup procedure as for nonpolar columns. For a column with nonaqueous mobile phases, use the following scheme:</p> <ul style="list-style-type: none"> Flush with 50ml chloroform Flush with 50ml methylene chloride Flush with 50ml 2-propanol Flush with 25ml methylene chloride Flush with 25ml mobile phase Evaluate column performance
<p>Nonpolar Column (C1, C4, C6, C8, C18, phenyl, or diphenyl column, SUPELCOSIL® LC-PAH or DB column)</p> <p>Water Soluble Samples</p> <ul style="list-style-type: none"> Flush with 50ml hot (40-60°C) distilled water Flush with 50ml methanol Flush with 50ml acetonitrile Flush with 30ml tetrahydrofuran Flush with 25ml methanol Flush with 25ml mobile phase Evaluate column performance according to conditions specified by the manufacturer 	<p>Ion Exchange Column (strong or weak anion or cation exchange, amino column)</p> <p>Most analyses involving ion exchange systems use ionic mobile phases. Compounds that may affect column performance are usually insoluble or only slightly soluble in water. The following procedure should be sufficient to remove these compounds.</p> <ul style="list-style-type: none"> Flush with 50ml hot (40-60°C) distilled water Flush with 50ml methanol Flush with 50ml acetonitrile Flush with 25ml methylene chloride Flush with 25ml methanol Flush with 25ml mobile phase Evaluate column performance
<p>Water Insoluble Samples</p> <ul style="list-style-type: none"> Flush with 50ml isopropanol Flush with 50ml tetrahydrofuran Flush with 50ml methylene chloride Flush with 50ml hexane Flush with 25ml isopropanol Flush with 25ml mobile phase Evaluate column performance 	<p>Columns Used for Protein/Biopolymer Analyses (C4, C8, C18, other RP column used for proteins, diphenyl column)</p> <ul style="list-style-type: none"> Flush with 100ml distilled water Flush with 50ml 0.1% trifluoroacetic acid Flush with 50ml isopropanol Flush with 50ml acetonitrile Flush with 50ml distilled water Flush with 50ml mobile phase Evaluate column performance